



# Bio-signal research

Guger Technologies OG (g.tec)

Julita de la Vega Arias



ACHI 2012 - January 30 - February 4, 2012 - Valencia, Spain



# 1. Guger Technologies OG (g.tec)

## Company fields


- bio-engineering, medical electronics
- developing and offering hard- and software products for biosignal research  
(**single cell activity, EEG, ECoG, ECG, EMG, EOG, ...**)
- user specific adaptations and developments
- mainly based on rapid prototyping environment under *MATLAB/Simulink*



# 1. Guger Technologies OG (g.tec)

## Company description

- private company
- inter-disciplinary team (biomedical- , telematics engineers, bio-medical engineers, software engineers, macrotronic engineers, etc)
- customers: universities, university hospitals, R&D departments, industry





# 1. Guger Technologies OG (g.tec)

## Our Offices

- Guger Technologies OG,  
Herbersteinstrasse 60, 8020 Graz, Austria
- g.tec medical engineering GmbH,  
Sierningstrasse 14, 4521 Schiedlberg, Austria
- g.tec medical engineering Spain,  
Calle Pallars 74-76, 3-1, Barcelona, Spain

## Contact

**phone:** +43 316 675106, **fax:** +43 316 675106-39

**Web:** [www.gtec.at](http://www.gtec.at), [office@gtec.at](mailto:office@gtec.at)






# 1. Guger Technologies OG (g.tec)

- World-wide supplier in more than 60 countries
- Most important markets: USA, Japan, Europe
- Currently several 100 institutions using the technology
  - Wadsworth Center, Albany, NY
  - Riken Institute, Japan together with Toyota, Japan
  - NASA, USA
  - MIT, Harvard, Johns Hopkins, Stanford, Cambridge...

## 2. Participation in EC projects

- **ReNaChip** Rehabilitation of a discrete sensory motor learning function
- **Sm4all** Brain-Computer Interface for smart home control
- **RGS** faster recovery from stroke with games
- **BrainAble** BCI with VR, smart homes and social networks
- **Decoder** Assessment/communication in non-responsive patients
- **CSI** Advance the state of the art in 3D medical-imaging platforms
- **BETTER** BCI for Stroke rehabilitation using rehabilitation robots
- **VERE** human body and surrogate representation in virtual and physical reality
- **ALIAS** cognitive assistance and social interaction by robots





### 3. Our co-operations partners

- **University of Barcelona, Spain**  
Mel Slater, Chris Groenegrass
- **IDIBAPS, Barcelona, Spain**  
Mavi Sanchez-Vives
- **University College London (UCL), UK**  
Anthony Steed, Angus Antely,
- **University of Technology Graz, Austria**  
Robert Leeb, Gert Pfurtscheller
- **Wadsworth Center, New York, USA**  
Gerwin Schalk, Eric Sellers
- **Tel Aviv University**  
Matti Mintz





## 4. Hardware Development

EEG, ECoG, EMG, EOG, ECG amplifiers 1-128 channels

- A. USB based biosignal amplifier
- B. Mobile device
- C. Stand-alone device





**We built an amplifier that can now be used not only for EEG, but also for implanted electrodes in humans, spike recordings and in vitro setups.**



**g-USBamp 3.0**  
USB BIO SIGNAL AMPLIFIER

explore the colors of science.

16 DC-channels per unit.  
4 independent grounds:  
record any type of signal:  
(EEG/ECG/ECG/EMG/EOG/...)  
connect various sensors.  
stack units for 32/48/64/... ch.

24-bit resolution with  
simultaneous sampling of all  
channels with up to 38.4 kHz.  
digital signal filtering and  
preprocessing.  
connect via USB 2.0

works with passive and with  
active electrodes.  
8 digital trigger inputs/unit.  
4 digital output/unit.  
new simplified  
synchronization of units.

internal digital bandpass and  
notch filters.  
built-in calibration unit and  
impedance checking.  
power supply and accumulator  
package available.

easy configuration and setup  
via the software.  
high-speed online data  
processing for MATLAB/  
SIMULINK and for LabVIEW  
available.

used by leading research  
groups for non-invasive and  
for invasive recording.  
available as part of a complete  
BCI (Brain-Computer Interface)  
research system.

CE-certified medical device.  
safety class: II  
conformity class: IIIa  
type of applied part: CF  
FDA approved device.

contact: g.tec medical engineering GmbH, Sierningstr. 14, 4521 Schiedberg, AUSTRIA, phone: +43 7251 22240, fax: +43 7251 22240 39, e-mail: office@gtec.at, web: www.gtec.at

[www.gtec.at](http://www.gtec.at)

# g. SAHARA

ACTIVE DRY ELECTRODE SYSTEM

active dry EEG electrode system



## **g.tec developed active and dry EEG electrode**

g.SAHARA electrode system (patent pending) consists of an 8 pin electrode made of a special golden alloy

Pins have sufficient length to reach through the hair to the skin

Golden alloy and the 8 pins reduce the electrode-skin impedance

Electrode itself can be connected with a clip to the active electrode system on top of it



**g·Hamp**  
MULTICHANNEL AMPLIFIER



80/144/256  
channels

**g·tec**  
GUGER  
TECHNOLOGIES

- 256 channels per unit with 4 multi-pole connectors at front side
- 1 ADC per channel for synchronization
- 24 Bit ADC per channel for perfect resolution
- Input range of +/-250mV to measure from EMG, ECG, ECoG, Spikes, EEG down to EPs with only one setting. Therefore no saturation of the amplifier possible
- DC – 6 kHz input frequency range
- CE approved
- C API, MATLAB API, Simulink driver Linux and Windows version
- OS 32 Bit and 64 Bit
- Different headboxes for passive or active
- Electrodes or headstages
- 16 digital inputs/outputs for triggering with external devices
- Impedance check included



[www.gtec.at](http://www.gtec.at)

## 5. Software Programming Environment

### A. C++ Application Program Interface (API)

- allows to integrate amplifiers into own software under Windows and Linux

### B. MATLAB API

- integrate amplifiers into MATLAB data acquisition and analysis programs
- access to all toolboxes (Signal Processing, Neural Networks,...)
- access to user written M-files

### C. Simulink Highspeed on-line Processing

- amplifier device driver block under Simulink
- copy the block into Simulink model and connect the signal processing (S-functions) and paradigm blocks (MATLAB code)
- just exchange the amplifier device driver and work with the same signal processing blocks

### D. LabView

- amplifier device driver block under LabView
- use standard LabView blocks for analysis



**THANKS!!!!**

